

CLEANING GUIDE FOR

HCL® Paralytic Box w/ 3 Compartments (#21330)

Materials:

Acrylic

Recommended Disinfectant Options:

19954 PREempt™ One-Step Surface Cleaner and Disinfectant Wipes, 6 x 7, Canister

19954-31 PREempt™ One-Step Surface Cleaner and Disinfectant Wipes, 6 x 7, Case

Alternate Disinfectant Options:

3% Hydrogen Peroxide based cleaners

15% Sodium Hypochlorite based cleaners

Follow the cleaner's instructions and adhere to the specific Standard Operating Procedure (SOP) established by your healthcare system.